

Essential Energy

Valleylab™ LS10

LS SERIES SINGLE CHANNEL VESSEL SEALING GENERATOR

3rd generation of LigaSure™ technology

INTELLIGENT TISSUEFFECT™ PLUS TECHNOLOGY¹

Applies appropriate energy to targeted tissue, for an effective seal with minimal thermal spread in vessels up to 7 mm and tissue bundle.

FAST SEALING

Seals in 2-4 seconds¹



INTUITIVE CONTROL PANEL

Single power on/off button, single LigaSure™ instrument plug-in port.³

SMART™ CONNECTOR

Recognizes which LigaSure™ instrument is being used and automatically configures energy settings for quick, consistent results.

COMPACT & LIGHTWEIGHT

Uses minimal space and is easily portable.³
Size: 30.0 x 10.5 x 37.7 cm.
Weight: 5.5 kg.

LigaSure™ technology is now accessible everywhere. Benefit from its simple, plug-and-play design as well as its portability when transporting the generator between operating rooms. The Valleylab™ LS10 generator delivers the next generation of LigaSure™ technology, allowing anyone in the Operating Room to reach beyond traditional ligation.